

PRODUCT FOCUS:

CHAIR OCCUPANCY SENSOR



TASK Community Care | Product Focus

Chair Occupancy Sensor



Part of the TASK CareLink Falls and Memory Support Packages, for service users who have a substantial or critical need. Please contact us for further details on this or other TASK CareLink Call Packages

For more information on the Chair Occupancy Sensor contact us at:

TASK Community Care
Unit 1, Block 3
CityNorth Business Park
Stamullen
County Meath
IRELAND

info@taskltd.ie
Telephone: +353 (0)1 8435889
www.taskltd.ie

TASK Community Care
The Innovation Centre
Northern Ireland Science Park
Queen's Road
BELFAST
BT3 9DT

info@taskltd.com
Telephone: 0845 304 5535
www.taskltd.com

Overview

The Chair Occupancy Sensor is placed out of sight under the cushion of a chair. It can send an alert immediately when an individual leaves the chair, indicating a possible fall. Alternatively, it can be programmed to send the alert only if the individual does not return within a reasonable time.

This device can be used on its own or installed as part of an integrated telecare system.

How it works

The Chair Occupancy Sensor can operate in a number of ways specific to an individual's unique set of requirements:

1. The device can automatically sound an attached buzzer as soon as a person leaves the chair. An onsite carer can then ensure the individual is safe.
2. Similarly, it can be used with a radio trigger which means the buzzer will sound up to 200meters away, in-house or to a neighbouring home.
3. Finally, with any of the above options or as an alternative, the system can be used with one of our CareLink two-way speech units. This means staff at our 24/7 CareLink Monitoring Centre can talk with the individual through the CareLink two-way speech unit and arrange for help if needed.

An additional option is to have the device automatically switch on a light when the person gets off the chair. This can help to minimise the risk of having a fall if the room is not adequately lit.

Alternatively we can supply fully compatible replacement Chair Occupancy Sensors for your existing NurseCall system.

Benefits

Reduces the possibility of an individual falling and being left lying on the floor indefinitely. Supports independent living.

Warranty

Chair Occupancy Pad – 1 Year Warranty
CareLink two-way speech unit – Lifetime Warranty